200 MHz Photoreceiver with Si-PIN Photodiode



The picture shows model HCA-S-200M-SI-FST

Features Si-PIN photodiode, 0.8 mm active diameter Bandwidth DC - 200 MHz Amplifier transimpedance gain 2.0 × 10⁴ V/A Max. conversion gain 1.1 × 10⁴ V/W @ 800 nm Spectral range 320 - 1000 nm Free-space input 1.035"-40 threaded, easily convertible to fiber optic input (FC and FSMA) with optionally available screw-on adapters Fiber optic input also available as permanently mounted FC-input UNC 8-32 and M4 tapped holes for mounting on standard posts with metric and imperial thread **Applications Spectroscopy** Fast pulse and transient measurements **Optical triggering** Optical front-end for oscilloscopes, A/D converters and HF lock-in amplifiers Block Diagram Rf OPTICAL **INPUT** Buffer VOLTAGE I/V amplifier OUTPUT Bias Offset

nulling

Intended Use

The HCA-S-200M-SI photoreceiver consists of an Si photodiode and a subsequent low-noise fixed gain transimpedance amplifier. It is designed for fast conversion of small optical signals into equivalent output voltages. Operation is mostly self-explanatory. If in doubt, consult this document or contact support@femto.de.

For safe operation, please refer to the damage thresholds specified in the "Absolute Maximum Ratings", "Temperature Range" and "Power Supply" sections of this document.

The operating environment must be free of smoke, dust, grease, oil, condensing moisture, and other contaminants that could affect the operation or performance.

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

200 MHz Photoreceiver with Si-PIN Photodiode

Available Versions

HCA-S-200M-SI-FST



1.035"-40 threaded flange with internally threaded coupler ring (outer diameter 30 mm) for free space applications. Compatible with many optical standard accessories and for use with various types of fiber connector adapters.

Optionally available:

Fiber adapters PRA-FC, PRA-FCA and PRA-FSMA. With the relative large 0.8 mm dia. photodiode installed in the HCA-S-200M-SI input coupling is not critical. However, standard SM 9/125 fibers (PC or APC) with low numerical aperture (NA) are recommended for ensuring near 100% coupling efficiency.

HCA-S-200M-SI-FC



Fix/permanent FC fiber connector for high coupling efficiency and excellent conversion gain accuracy.

Related Models

HCA-S-200M-IN-FST

HCA-S-200M-IN-FC

InGaAs-PIN, Ø 0.3 mm, 900 - 1700 nm free space input, 1.035"-40 threaded flange

InGaAs-PIN, integrated ball lens, 900 - 1700 nm FC fiber connector (fix/permanent)

Available Accessories

PRA-FCA PRA-FSMA







Fiber-adapter with external 1.035"-40 thread

PRA-PAP



Alternative mounting option: Post adapter plate, easy to mount on FEMTO photoreceiver series OE, FWPR, PWPR, HCA-S and LCA-S

PS-15-25-L



Power Supply Input: 100 – 240 VAC Output: ±15 VDC

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

200 MHz Photoreceiver with Si-PIN Photodiode

 $V_S = \pm 15 \text{ V}, T_A = 25 \,^{\circ}\text{C}, \text{ output load impedance } 50 \,\Omega,$ Specifications Test conditions

warm-up 20 minutes (min. 10 minutes recommended)

Gain Transimpedance gain 2.0×10^4 V/A (@ output load 50 Ω)

> ±1 % (electrical) Gain accuracy

 1.1×10^4 V/W typ. (@ 800 nm, output load 50 Ω) Conversion gain

Frequency Response Lower cut-off frequency

Upper cut-off frequency (-3 dB) 200 MHz (±10 %)

Gain flatness $\pm 1 dB$

Time Response Rise/fall time (10 % - 90 %) 1.8 ns

Input Noise equivalent power (NEP) 9.4 pW/√Hz (@ 800 nm, 10 MHz)

Optical saturation power 110 µW (for linear amplification, @ 800 nm) Input offset compensation range ±100 μA, adjustable by offset potentiometer

Detector Si-PIN photodiode Detector

Active area Ø 0.8 mm Spectral range 320 - 1000 nm

0.55 A/W typ. (@ 800 nm) Max. sensitivity

Output Output voltage range $\pm 1.2 \text{ V}$ (@ 50 Ω output load)

for linear operation and low harmonic distortion

Max. output voltage range $\pm 1.7 \text{ V } (@ 50 \Omega \text{ load})$

Output impedance 50 Ω (terminate with 50 Ω load)

Output noise 3 mV RMS (20 mV peak-peak) typ. (@ 50 Ω load,

no signal on detector, measurement bandwidth 500 MHz)

Optical Input Connector 1.4305 stainless steel, nickel-plated Material FST flange

Material FST coupler ring 1.4305 stainless steel, glass bead blasted

Material FC receptacle nickel silver

Power Supply Supply voltage ±15 V (±14.5 V ... ±16.5 V)

> Supply current ±50 mA (depends on operating conditions,

recommended power supply capability min. ±150 mA)

Case Weight 209 g (0.46 lbs) HCA-S-200M-SI-FST incl. coupler ring

188 g (0.41 lbs) HCA-S-200M-SI-FC

Material AlMg4.5Mn, nickel-plated

Temperature Range Storage temperature -30 °C ... +85 °C

Operating temperature 0 °C ... +60 °C

Absolute Maximum Ratings Optical input power (CW) 20 mW

Power supply voltage ±20 V

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

200 MHz Photoreceiver with Si-PIN Photodiode

Connectors

Input

HCA-S-200M-SI-FST

1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories


HCA-S-200M-SI-FC

FC fiber optic connector

HCA-S-200M-SI-FC FC fiber optic connector (fix/permanent, FC/PC and

Output BNC jack (female)

Power supply LEMO® series 1S, 3-pin fixed socket (mating plug type: FFA.1S.303.CLAC52)



Scope of Delivery HCA-S-200M-SI, internally threaded coupler ring (FST version only), LEMO® 3-pin connector, datasheet, transport package


Ordering Information HCA-S-200M-SI-FST 1.035"-40 threaded flange for free space applications and

for use with various types of optical standard accessories.

HCA-S-200M-SI-FC FC fiber optic connector

(fix/permanent, FC/PC and FC/APC compatible).


Spectral Response



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

200 MHz Photoreceiver with Si-PIN Photodiode

Typical Performance Characteristics



PD-HCA-S-200M-Si-bw_R01

Noise spectrum

* RBW 1 MHz

Att 0 dB * VBW 1 kHz Noise2 292.328379 nV/√Hz Ref 22.4 mV SWT 800ms 254.400000000 MHz 103.266894 nV/√Hz Noise1 10 mV 10.00000000 MHz 1 mV М2 140.uv 10 μV 1 μV-

Start 0.0 Hz

Stop 400.0 MHz

PD-HCA-S-200M-Si-noise-R01


Note: spectral noise data is measured at the amplifier output with no signal on the photodiode. To determine the spectral input noise divide the measured output noise by the amplifier conversion gain. Conversion gain (V/W) = amplifier gain (V/A) × photo sensitivity (A/W).

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

200 MHz Photoreceiver with Si-PIN Photodiode


Typical Performance Characteristics (continued)

Pulse response to square wave input signal (with 16 times averaging)



PD-HCA-S-200M-Si pulse-2ns_R01


Large signal response output signal for 100 MHz, 100 μW modulated optical input signal (with 4 times averaging)




SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

200 MHz Photoreceiver with Si-PIN Photodiode

Typical Performance Characteristics (continued) Small signal response output signal for 1.5 μW modulated optical input signal, 1 MHz square wave, without averaging



Small signal response output signal for 1.5 μW modulated optical input signal, 1 MHz square wave, with 64 times averaging





SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

200 MHz Photoreceiver with Si-PIN Photodiode

Dimensions

HCA-S-200M-SI-FST (1.035"-40 threaded free space input)





DZ-HCA-S_FST_R1

all dimensions in mm unless otherwise noted

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY


Datasheet


HCA-S-200M-SI

200 MHz Photoreceiver with Si-PIN Photodiode

Dimensions (continued)

HCA-S-200M-SI-FC (FC fiber optic connector)





DZ-HCA-S_FC_R1

all dimensions in mm unless otherwise noted

FEMTO Messtechnik GmbH Klosterstr. 64 10179 Berlin · Germany Phone: +49 30 280 4711-0 Fax: +49 30 280 4711-11 Email: info@femto.de www.femto.de Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

© by FEMTO Messtechnik GmbH · Printed in Germany

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E T O