


200 MHz Photoreceiver  
with InGaAs-PIN Photodiode



The picture shows model HCA-S-200M-IN-FST

Features	<ul style="list-style-type: none"><li>• InGaAs-PIN photodiode</li><li>• Bandwidth DC – 200 MHz</li><li>• Amplifier transimpedance gain <math>2.0 \times 10^4</math> V/A</li><li>• Max. conversion gain <math>1.9 \times 10^4</math> V/W @ 1550 nm</li><li>• Spectral range 900 – 1700 nm</li><li>• Free-space input 1.035"-40 threaded</li><li>• Fiber optic input available as permanently mounted FC-input</li><li>• UNC 8-32 and M4 tapped holes for mounting on standard posts with metric and imperial thread</li></ul>
Applications	<ul style="list-style-type: none"><li>• Spectroscopy</li><li>• Fast pulse and transient measurements</li><li>• Optical triggering</li><li>• Optical front-end for oscilloscopes, A/D converters and HF lock-in amplifiers</li></ul>
Block Diagram	<p>The block diagram illustrates the internal circuitry of the photoreceiver. It begins with an "OPTICAL INPUT" which is coupled to an InGaAs photodiode. The photodiode is biased by a "Bias" voltage. The current generated by the photodiode flows into the inverting input of a transimpedance amplifier, represented by a triangle labeled "I/V". A feedback resistor, labeled "Rf", connects the output of the I/V amplifier back to its inverting input. An "Offset nulling" input is also shown connected to the non-inverting input of the I/V amplifier. The output of the I/V amplifier is then fed into a "Buffer amplifier", which finally produces the "VOLTAGE OUTPUT".</p> <p>BS01-HCA-S_R02</p>
Intended Use	<p>The HCA-S-200M-IN photoreceiver consists of an InGaAs photodiode and a subsequent low-noise fixed gain transimpedance amplifier. It is designed for fast conversion of small optical signals into equivalent output voltages. Operation is mostly self-explanatory. If in doubt, consult this document or contact <a href="mailto:support@femto.de">support@femto.de</a>.</p> <p>For safe operation, please refer to the damage thresholds specified in the "Absolute Maximum Ratings", "Temperature Range" and "Power Supply" sections of this document.</p> <p>The operating environment must be free of smoke, dust, grease, oil, condensing moisture, and other contaminants that could affect the operation or performance.</p>

## 200 MHz Photoreceiver with InGaAs-PIN Photodiode

### Available Versions

#### HCA-S-200M-IN-FST



1.035"-40 threaded flange with internally threaded coupler ring (outer diameter 30 mm) for free space applications, compatible with many optical standard accessories

#### HCA-S-200M-IN-FC



Fix/permanent FC fiber connector for high coupling efficiency and excellent conversion gain accuracy

### Related Models

#### HCA-S-200M-SI-FST

Si-PIN,  $\varnothing$  0.8 mm, 320 – 1000 nm  
free space input, 1.035"-40 threaded flange

#### HCA-S-200M-SI-FC

Si-PIN,  $\varnothing$  0.8 mm, 320 – 1000 nm  
FC fiber connector (fix/permanent)

### Available Accessories

#### PRA-PAP



Alternative mounting option:  
Post adapter plate, easy to mount  
on FEMTO photoreceiver series OE,  
FWPR, PWPR, HCA-S and LCA-S.

#### PS-15-25-L



Power Supply  
Input: 100 – 240 VAC  
Output:  $\pm 15$  VDC

### Specifications

#### Test conditions

$V_S = \pm 15$  V,  $T_A = 25$  °C, output load impedance 50  $\Omega$ ,  
warm-up 20 minutes (min. 10 minutes recommended)

#### Gain

Transimpedance gain  
Gain accuracy  
Conversion gain

$2.0 \times 10^4$  V/A (@ output load 50  $\Omega$ )  
 $\pm 1$  % (electrical)  
 $1.9 \times 10^4$  V/W typ. (@ 1550 nm, output load 50  $\Omega$ )

#### Frequency Response

Lower cut-off frequency  
Upper cut-off frequency (–3 dB)  
Gain flatness

DC  
200 MHz ( $\pm 15$  %)  
 $\pm 1$  dB

#### Time Response

Rise/fall time (10 % – 90 %)

1.8 ns

#### Input

Noise equivalent power (NEP)  
Optical saturation power  
Input offset compensation range

5.2 pW/ $\sqrt{\text{Hz}}$  (@ 1550 nm, 10 MHz)  
60  $\mu\text{W}$  (for linear amplification, @ 1550 nm)  
 $\pm 100$   $\mu\text{A}$ , adjustable by offset potentiometer


## 200 MHz Photoreceiver with InGaAs-PIN Photodiode

### Specifications (continued)

Detector	Detector	InGaAs-PIN photodiode
	Active area (FST version)	Ø 0.3 mm
	Active area (FC version)	integrated ball lens
		suitable for fibers up to 62.5 µm core diameter
	Spectral range	900 – 1700 nm
	Max. sensitivity	0.95 A/W typ. (@ 1550 nm)
Output	Output voltage range	±1.2 V (@ 50 Ω output load) for linear operation and low harmonic distortion
	Max. output voltage range	±1.7 V (@ 50 Ω output load)
	Output impedance	50 Ω (terminate with 50 Ω load)
	Output noise	4.5 mV RMS (30 mV peak-peak) typ. (@ 50 Ω load, no signal on detector, measurement bandwidth 500 MHz)
Optical Input Connector	Material FST flange	1.4305 stainless steel, nickel-plated
	Material FST coupler ring	1.4305 stainless steel, glass bead blasted
	Material FC receptacle	nickel silver
Power Supply	Supply voltage	±15 V (±14.5 V ... ±16.5 V)
	Supply current	±60 mA (depends on operating conditions, recommended power supply capability min. ±150 mA)
Case	Weight	209 g (0.46 lbs) HCA-S-200M-IN-FST incl. coupler ring 188 g (0.41 lbs) HCA-S-200M-IN-FC
	Material	AlMg4.5Mn, nickel-plated
Temperature Range	Storage temperature	–30 °C ... +85 °C
	Operating temperature	0 °C ... +60 °C

Absolute Maximum Ratings	Optical input power (CW)	10 mW
	Power supply voltage	±20 V

Connectors	Input	HCA-S-200M-IN-FST	1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories
		HCA-S-200M-IN-FC	FC fiber optic connector (fix/permanent, FC/PC and FC/APC compatible)
	Output	BNC jack (female)	
	Power supply	LEMO® series 1S, 3-pin fixed socket (mating plug type: FFA.1S.303.CLAC52)	





Pin 1: +15 V  
Pin 2: –15 V  
Pin 3: GND

Scope of Delivery	HCA-S-200M-IN, internally threaded coupler ring (FST version only), LEMO® 3-pin connector, datasheet, transport package
-------------------	---

200 MHz Photoreceiver  
with InGaAs-PIN Photodiode


Ordering Information	HCA-S-200M-IN-FST	1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories.
	HCA-S-200M-IN-FC	FC fiber optic connector (fix/permanent, FC/PC and FC/APC compatible).



200 MHz Photoreceiver  
with InGaAs-PIN Photodiode

Typical Performance  
Characteristics (continued)

Noise spectrum




PD-HCA-S-200M-IN-noise\_R01

Note: spectral noise data is measured at the amplifier output with no signal on the photodiode.  
To determine the spectral input noise divide the measured output noise by the amplifier  
conversion gain.

Conversion gain (V/W) = amplifier gain (V/A) × photo sensitivity (A/W).

Marker	frequency	output noise	resulting input noise (NEP)
1	10 MHz	93 nV/√Hz	4.9 pW/√Hz (@ 1550 nm)

Pulse response to square wave input signal  
(with 16 times averaging)




PD-HCA-S-200M-IN-pulse-2ns\_R01


200 MHz Photoreceiver  
with InGaAs-PIN Photodiode

Typical Performance  
Characteristics (continued)

Large signal response  
output signal for 100 MHz, 55  $\mu$ W modulated optical input signal  
(with 4 times averaging)




Small signal response  
output signal for 1.2  $\mu$ W modulated optical input signal, 1 MHz square wave, without averaging



200 MHz Photoreceiver  
with InGaAs-PIN Photodiode


Typical Performance  
Characteristics (continued)

Small signal response  
output signal for 1.2  $\mu$ W modulated optical input signal, 1 MHz square wave,  
with 64 times averaging



Dimensions

HCA-S-200M-IN-FST (1.035"-40 threaded free space input)



## Dimensions (continued)

Technical drawing showing the 1000 Series Motor Mounting Dimensions. The drawing includes three views: a top view, a side view, and a front view, along with a detail of the mounting hole.

**Top View Dimensions:**

- Total width: 94
- Width to center of terminal block: 87
- Width to center of mounting hole: 74

**Side View Dimensions:**

- Height: 51

**Front View Dimensions:**

- Total height: 27.3
- Height to center of terminal block: 14.7

**Detail Dimensions:**

- Mounting hole diameter:  $\varnothing 3.2$

D7-HCA-S EC: B1

all dimensions in mm unless otherwise noted

Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

**F E M ▶ T O**